




From: Separation Science e-Learning <noreply@sepscience.com>
Sent: Monday, December 03, 2012 1:07 PM
To: Hanchett, James (DPH)
Subject: Today in Separation Science - Latest Issue of Separation Science



Web Version Forward Unsubscribe



Volume 4 Issue 13



Read The Latest Issue of Separation Science, which includes articles on:

GPC/SEC Troubleshooting #7 Elution Volume/Retention Time Shift
Q: I have measured my calibration standards and samples and the elution volume/retention time has significantly changed. All signals are shifted to higher elution volume. What could be the reason for that?
[Click for PDF>>](#) [View Archive>>](#)

A Generic Method Development Strategy for Liquid Chromatography
We propose a novel method development (MD) strategy that still works when the compounds are not available individually, when no mass spectrometer is available and/or when the number of compounds is not known. Since MD usually requires a lot of trial and error it can take up several days or even weeks work, computer-assisted MD can, therefore, be very helpful. The existing MD strategies are either search-based (e.g. Simplex) or model-based (e.g. Drylab, Chromsword), or both (Design of Experiments). Model-based strategies require the modelling of the retention properties of all sample compounds and, therefore, require peak tracking (using DAD- or MS-spectra or injecting each compound individually). Moreover, modelling errors can lead to wrong solutions.
[Click for PDF>>](#) [View Archive>>](#)

Analysis of Multiresidue Pesticides from Food Using the QuEChERS Sample Preparation Approach, LC-MS-MS and GC-MS Analysis
Analysis of food products is challenging due to the variety and complexity of both the matrices and the compounds of interest. Sample preparation and downstream analysis require careful consideration to ensure method robustness as well as accurate and precise quantification. In this study we explore the analysis of multiple pesticide residues in spinach samples. This article will walk through the step-by-step process of developing the analytical method, from sample preparation to analysis, best suited to the data requirements.
[Click for PDF>>](#) [View Archive>>](#)

Multidimensional Gas Chromatography — A Technology For Today's Analytical Challenges
Today, the discerning chromatographer should expect that with modern instrumentation, the gas chromatography (GC) method should deliver an advanced capability to assess the chemical constituents of a sample, with precision, and a significant coverage of sample components present in trace to major amounts. Various strategies are available to meet this goal, albeit with varying degrees of success, and this can be gleaned by dissecting the different components of the GC method — the injector, the separating column and the detector.
[Click for PDF>>](#) [View Archive>>](#)


Time of Flight Mass Spectrometry for The Modern Lab: Qualitative Analysis, Quantitative Measurement and Rapid Screening and Identification
The latest Time-of-Flight Mass Spectrometers (TOF-MS) offer scientists a platform for quantification combined with the power of exact mass. Laboratories use of MS technology has expanded as the needs for fast accurate results have increased. These challenges range from measuring purity of drugs, identification of contaminants in food, to screening and confirmation of illicit drugs. This article explains the use of a TOF-MS system that provides an MS technology that enables the scientist to rapidly screen, identify and quantify compounds with accurate mass, screening speed and ease of use.
[Click for PDF>>](#) [View Archive>>](#)

FEATURED APPLICATIONS

Optimizing performance of the Agilent 1290 Infinity LC System using 1-mm id Columns
Agilent Technologies
[Click to read>>](#)

Extraction of Synthetic Cannabinoids (SPICE) and Metabolites from Urine, Plasma and Whole Blood using ISOLUTE SLE+ Prior to LC-MS/MS Analysis
Biotage
[Click to read>>](#)

Unattended Automated Standard Addition and Headspace Analysis for the Quantitative Determination of VOCs in Food Packaging
Thermo Scientific



ON-DEMAND WEBINARS

Why Inertness Matters In Gas Phase Analyses
Presenter: Ken Lynam
[Click here>>](#)

Trouble-free Food Testing: Easy Technologies to Improve Uptime and Productivity
Presenters: Dr Jason Link & Mike Chang
[Click here>>](#)

Trends in Method Development for HPLC and UHPLC - the Movement Toward Using Solid-Core Particles
Presenter: Richard A. Henry
[Click here>>](#)

How to get Rugged Results in Pesticide Analysis using Triple Quadrupole GC-MS/MS?
Presenter: Katerina Mastovska
[Click here>>](#)

Making the Undetectable Detectable — Derivatization in GC Analysis
Presenter: Ruth Lewis
[Click here>>](#)

Pesticide Screener: A Comprehensive Approach to Multi-Target Pesticide Screening Using LC- Time-of-Flight Mass Spectrometry
Presenter: George McLeod
[Click here>>](#)

Using Variations in Solid Core Particle Diameter and Pore Size to Improve UHPLC and HPLC Separations
Presenter: Stephen Luke
[Click here>>](#)

Stationary Phase and Mobile Phase Selection for Hydrophilic Interaction Liquid Chromatography
Presenter: Dafydd Milton
[Click here>>](#)

Addressing the Needs of Bio-analysis with New SPE Technologies — SOLA
Presenter: Mike Oliver
[Click here>>](#)

Emerging Technologies, Low- & No-Prep Methods for Rapid MS Analyses in Juices
Presenter: Courtney Phillips
[Click here>>](#)

Fast LC Method Development Using Superficially Porous Particle Columns
Presenter: John Dolan
[Click here>>](#)

1

Folk_OIG_PRR_093947

[Click to read>>](#)

Online Analysis of Trace Anions in Power Plant Water
Matrices
Metrohm

[Click to read>>](#)

Molecular Weight Determination of LMWH SEC/MALS vs.
SEC/UV-RI
Wyatt Technology

[Click to read>>](#)

Dramatically Improve Existing 5 µm and 3 µm Fully
Porous Methods with Kinetex 5 µm Core-shell Technology
Phenomenex

[Click to read>>](#)

Automated USP and EP GPC/SEC Analysis of Low
Molecular Weight Heparin (LMWH)
PSS

[Click to read>>](#)

Fast Analysis of Water Pollutant Pesticides
Knauer

[Click to read>>](#)

Development of an SPME-GC-MS Procedure for the
Determination of Phthalate Esters in Ramen Noodle Flavor
Packets
Supelco/Sigma-Aldrich

[Click to read>>](#)

Published by Eclipse Business Media Ltd
Frederick House | Princes Court | Beam Heath Way | Nantwich | Cheshire CW5 6PQ | United Kingdom
20 Maxwell Road | #09-17 Maxwell House | Singapore 069113

Copyright © 2012 Eclipse Business Media Ltd. All rights reserved.

To sponsor an advert in this e-newsletter, contact david.grimes@sepscience.com
For editorial contributions, contact david.hill@supelco.com

This message was sent from Separation Science e-Learning to james.hanchett@state.ma.us. It was sent from: Eclipse Business Media Ltd, Frederick House, Princes Court, Beam Heath Way, Nantwich, Cheshire CW5 6PQ, United Kingdom. You can modify/update your subscription via the link below.



[Unsubscribe](#)